US

EN: Indication for use

Medium for culture of embryos from fertilisation to the blastocyst stage and for embryo transfer.

Contraindications

Gx-TL™ contains gentamicin and acetylcysteine. Do not use in patients with known hypersensitivity/allergy to any of the components.

Product Description

SUPPLEMENTED WITH HSA.

Gx-TL is a bicarbonate buffered medium containing human serum albumin, hyaluronan, acetylcysteine, and gentamicin as an antibacterial agent.

Ready to use after equilibration at +37 °C and 6 % CO2 atmosphere at sea level.

The summary of safety and clinical performance can be found at [www.vitrolife.com](#)

Storage instructions and stability

Store dark at +2 to +8 °C

Gx-TL is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the container. Opening of the bottle shall be limited to the necessary amount to remove the necessary amount of medium. Discard excess media no later than two weeks after first opening.

Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes fertility treatment.

Record the opening date on the bottle.

Gx-TL is designed for culture in a CO2 incubator.

Gx-TL shall be equilibrated at 6 % CO2 and +37 °C for no less than 6 hours and until correct pH has been attained before use. The dish shall not be kept in incubator for more than 7 days.

To obtain desired pH, the incubator CO2 level may have to be adjusted taking into account incubator characteristics and environmental conditions such as altitude.

To minimize evaporation, immediately cover droplets of Gx-TL with oil for ART procedures. Cuitluring without oil requires humidified environment.

Culture

Prepare dishes with droplets of Gx-TL (≥ 25 µL) under oil, or use a large volume (≥ 500 µL) with or without oil cover. For group culture, a minimum of 10 µL Gx-TL per embryo shall be used, to minimize risk of depletion of nutrients. Equilibrate according to above.

For time-lapse medium, prepare dishes with Gx-TL medium volumes as recommended by the manufacturer.

Following two washing steps in Gx-TL, transfer the fertilized oocytes or embryos to the Gx-TL culture medium and return the dish to the incubator immediately.

Embryo transfer

Prepare and equilibrate a dish with Gx-TL and perform embryo transfer according to laboratory procedures. Rinsing of the embryo transfer catheter with Gx-TL before aspiration of one or more embryos is recommended. Following the embryo transfer procedure, make a final microscopic examination of the catheter.

Specifications

Aseptically filled

Mouse Embryo Assay, 1-cell [% embryos developed to expanded blastocyst at 96 hours] ≥ 80

Bacterial endotoxins (LAL assay) [EU/mL] < 0.25

Osmaaliol tested

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use Gx-TL if it appears cloudy.

Gx-TL contains human serum albumin. Caution: All blood products should be treated as potentially infectious. Source material from this product was derived was found negative when tested for antibodies to HIV, HBC, HCV, and HTLV III and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

To avoid contamination VitroLife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including VitroLife's IVF media, have not been determined and are uncertain.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Symbols glossary

*Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021)

Symbol	Title	Reference and symbol ID
	Consult instructions for use	*Ref no: 5.4.3
	Caution	*Ref no: 5.4.4
	Manu-facturer	*Ref no: 5.1.1
	Date of manufacture	*Ref no: 5.1.3
	Use-by date	*Ref no: 5.1.4
	Batch code	*Ref no: 5.1.5
	Catalogue	*Ref no: 5.1.6
	Sterilized using aseptic processing	*Ref no: 5.2.2
	Single sterile barrier	*Ref no: 5.2.1.1
	Keep away from sunlight (light sources)	*Ref no: 5.3.2.2
	Temperature limit	*Ref no: 5.3.7
	Contains human blood or plasma derivatives	*Ref no: 5.4.6
	Contains a medicinal substance	*Ref no: 5.4.7
	Medical Device	*Ref no: 5.7.7
	CE Mark (Conformité Européen)	

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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